



# Release Notes – Version 3.0.0.8

---

*MODALIZER SDK - DICOM Toolkit*

*07 June 2020*

## Summary

number	Created on	Summary
B49	2018-09-05	Fixed Code crashing in RZDCX x64 with following code (Dyn) (3.0.0.4)
B99	2019-09-26	Fixed DCXELM. Value now supports assignment of arrays and performs more checks
B103	2019-11-14	Fixed Crash when setting CertSubjectName of server to a certificate in the personal store with no private key
B106	2019-12-22	Fixed DCXELM. Value of empty sequence throws unknown error
B110	2020-01-05	Fixed Clip Image wrong input frame number handling (3.0.0.7)
B113	2020-01-09	Fixed functionality of Jpeg2000 decoding in DCMOBJ (3.0.0.7)
B118	2020-02-19	Fixed MPPS Create may fail with corrupted SOP Instance UID in the request
B120	2020-02-23	Fixed Convert video to DICOM input file size limited to 2GB, now supported up to 4 GB
B122	2020-04-28	Fixed Assigning string value to sequence element causes access violation
B123	2020-05-06	Fixed DCXREQ. MoveAndStore() deadlocks



## Detailed Release Notes

### **B106- DCXELM. Value of empty sequence throws unknown error**

DCXELM with Value Representation SQ may only be assigned with DCXELMIterator