

Release Notes – Version 2.0.2.9

<u>#</u>	<u>Version Found</u>	<u>Created</u>	<u>Subsys</u>	<u>Title</u>	<u>Description</u>
353	2.0.2.8	Jan 24	DSRSVC	Correct Sequence Matching Behavior for Universal Sequence Matching	Sequence matching of empty sequence elements or single empty item should return the complete managed object of the sequence
354	2.0.2.8	Feb 05	RZDCX	Remove log message for tag not found or make it debug level message	This message should not be in Error level because it already throws an exception.
358	2.0.2.8	Feb 10	RZDCX	API for flushing log and writing log messages	Add API calls to DCXAPP: *: WriteToLog *: FlushLog
359	2.0.2.8	Feb 11	DSRSVC	Complete attributes of MWL to support GE LightSpeed CT	GE's LightSpeed CT has a strict attribute list for MWL that without all of them supported the device will not show worklist query results. The attribute list is described in the conformance statement: http://www3.gehealthcare.com/en/Products/Interoperability/DICOM/~/_media/Downloads/us/Product/Interoperability/DICOM/Computed-Tomography/GEHC-DICOM-Conformance_LightSpeed_DOC1246512_Rev1.pdf In Page 88 Fields considered Type 1 by the Modality Worklist Server include: *: (0010,0010), Patient Name *: (0010,0020), Patient ID *: (0020,000D), Study Instance UID *: (0040,0001), Scheduled Station AE Title *: (0040,0002), Scheduled Procedure Step Start Date 1 *: (0040,0003), Scheduled Procedure Step Start Time 1 *: (0040,0009), Scheduled Procedure Step ID *: (0040,1001), Requested Procedure ID 1 Start Date must be of the form YYYYMMDD, exactly eight numeric characters, and Start Time must be of the form HHMMSS, exactly six numeric

19 February 2013

					characters. Fields considered Type 2 by Modality Worklist Server include: *: (0008,0050), Accession Number *: (0008,0060), Modality *: (0008,0090), Referring Physician Name *: (0010,0030), Patient Date of Birth *: (0010,0040), Patient Sex *: (0010,1030), Patient Weight in kg *: (0010,2000), Medical Alerts *: (0010,2110), Contrast Allergies *: (0010,21C0), Pregnancy Status *: (0032,1032), Requesting Physician *: (0032,1070), Requested Contrast Agent *: (0038,0010), Admission ID *: (0038,0050), Special Needs *: (0038,0300), Current Patient Location *: (0038,0500), Patient State *: (0040,0006), Performing Physician *: (0040,0010), Scheduled Station Name *: (0040,0011), Scheduled Procedure Step Location *: (0040,0012), Pre-order Medication *: (0040,1003), Requested Procedure Priority *: (0040,1004), Patient Transport Arrangements *: (0040,3001), Confidentiality Constraint
--	--	--	--	--	--